


|  THE GOODELL LABORATORY | | |
|--|---|--------------|
| Author | Nathan Boles | Feb., 2009 |
| Title | Single Cell Real-Time PCR | |
| Introduction | This protocol describes single cell Real-time PCR. | |
| Materials | <ol style="list-style-type: none"> 1. 10 mM dNTP 2. 500ug/mL Random primer mix 3. dH₂O 4. 5X 1st strand buffer 5. Rnase inhibitor 6. NP40 7. 2X Taqman Master Mix 8. 18s Taqman probe 9. Taqman probes for your gene of interest (GOI) | |
| Protocol | | <i>Notes</i> |
| 1. | Prepare the stock random primer mix: <ol style="list-style-type: none"> 1. 40 uL 10mM dNTP 2. 20 uL 500ug/mL Random primer mix 3. 20 uL dH₂O Dilute 1:24 to make stock random primer mix | |
| 2. | Prepare lysis solution: <ol style="list-style-type: none"> 1. 608 uL dH₂O 2. 160 uL 5X 1st strand buffer 3. 16 uL of Rnase inhibitor 4. 4 uL NP40 5. 12 uL of stock random primer mix | |
| 3. | Pipet 8 uL of Lysis solution into each well of a 96 well plate | |
| 4 | Sort a single cell into each well, cover with optical cover sheet and bring back to lab | ... |

| | | |
|----|--|--|
| 5 | After sort, take 10 uL of Superscript II and mix with 190 uL of dH ₂ O. Now pipette 2 uL of mix into each well. Then do a quick spin of the plate to collect all liquid at the bottom of the well. | |
| 6 | Run plate on a PCR machine using a standard RT-PCR protocol | |
| 7 | <p>While RT-PCR is running prepare following master mix (does 20 wells):</p> <ol style="list-style-type: none"> 1. 275 uL 2X taqman master mix 2. 27.5 uL 18s taqman probe 3. 27.5 uL GOI probe | <i>20 individual cells per gene is reasonable, you can do more if you like just scale up the master mix.</i> |
| 8. | After RT-PCR reaction is done add 15 uL of master mix to each well | |
| 9. | Run real-time on ABI real-time system using standard real-time template | |