

	<h1 style="margin: 0;">THE GOODELL LABORATORY</h1>	
<p style="text-align: center;"><b>Author</b></p>	<p>Karen Lin</p>	<p>July 30, 2005</p>
<p style="text-align: center;"><b>Title</b></p>	<p><b>Titering MSCV Viruses with 3T3 Cells</b></p>	
<p style="text-align: center;"><b>Introduction</b></p>	<p>This protocol is to titer an ecotropic viral vector expressing GFP with FACS analysis of GFP expression. Target cell-lines and read-out of titering may be varied when distinct packaging method and viral construct are used.</p>	
<p style="text-align: center;"><b>Materials</b></p>	<p><b><i>Transduction Medium for 3T3 cells:</i></b></p> <p>DMEM (Invitrogen, 11965-092)          10% Calf Serum (Colorado Serum Co., Cat# 1334)          1X Pen/Strep          4 ug/ml polybrene (=Hexadimethrine bromide, Sigma, H9268)</p> <p><b><i>Culture Medium for 3T3 cells:</i></b></p> <p>DMEM (Invitrogen, 11965-092)          10% Calf Serum          1X Pen/Strep</p> <p><b>24-Well Tissue Culture Plates</b></p>	
<p style="text-align: center;"><b>Protocol</b></p>		<p><i>Notes</i></p>
<p style="text-align: center;"><b>1.</b></p>	<p><b>Day0</b>          (1) Plate <math>5 \times 10^4</math> 3T3 cells/well in a 24 well plate with 1 ml culture(C) medium.          (2) Incubate 37°C overnight.</p>	<p>keep healthy cell stock. Because cells will be under a huge stress from virus infection, cells for titering should be of healthy condition when loading viruses.</p>
<p style="text-align: center;"><b>2.</b></p>	<p><b>Day1</b>          (1) Replace medium with 950ul transduction (T) medium          (2) Thaw virus rapidly in 37°C waterbath.          (3) Make 10 fold dilution series of each virus, and apply 50ul of diluted sample onto cell culture.</p> <p style="text-align: center;"><math>10^0</math> = non-diluted thawed supernatant</p>	<p>Viral particles should be freshly thawed right before the infection.</p>

	$10^{-1}$ = 6 ul of non-diluted thawed supernatant + 54 ul T medium $10^{-2}$ = 6 ul of $10^{-1}$ + 54 ul T medium $10^{-3}$ = 6 ul of $10^{-2}$ + 54 ul T medium $10^{-4}$ = 6 ul of $10^{-3}$ + 54 ul T medium Leave one well for no virus control	
3.	<b>Day2 (24 hours after transduction)</b> Replace media with culture medium	
4.	<b>Day3 (48 hours after infection)</b> Trypsinize 3T3 cells, FACS analysis for GFP expression to evaluate titer For a well of ~1-10% expression in GFP <b>Titer = (%) GFP x (#) cell plated x dilution factor / 0.05 ml = viral particles/ml</b>	Typically with a decent virus production, we find the 1-10% eGFP bracket at the $10^{-2}$ dilution.

**Note:**

1. A freeze-thaw cycle would decrease the virus titer by ~10 fold. Therefore, it is very important to set an aliquot of viral supernatant for virus titering when harvesting viral particles out of 293T cells (the production of retrovirus).
2. Typically, for retroviruses produced from a 8kb vector, the titer is about  $5 \times 10^6$  virus/ml. For a virus from a 7kb vector, the titer can reach  $1 \times 10^7$  virus/ml.
3. For a good transduction of hematopoietic stem/progenitor cells, it is recommended to use virus titer no less than  $2 \times 10^6$  virus/ml. It is empirical that we found a larger volume of viral supernatant (exceeding a quarter of the total transduction volume) will decrease the transduction efficiency. For example, with a transduction of  $5 \times 10^5$  cells (in 1ml transduction medium), it is highly recommended to have a virus titer of no less than  $2 \times 10^6$  virus/ml.
4. If the titering results are not satisfying, it usually means the transfection reactions for virus production is not good. I would go back trouble-shooting the transfection efficiency.